

Certificate of Analysis

No. CA06165


Product: Eichrom's Sr Resin Cartridges

Lot #: SRSR6B

Eichrom's Sr Resin, Lot # SRSR6B has been manufactured according to Eichrom standard operating procedures. The lot of resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	133
Gravimetric Recovery	$\leq 105\%$	102 %
Chemical Recovery, Corrected	90-110 %	95.0 %



Jill Bryant

Quality Control Chemist

Prepared: December 20, 2006