

## Certificate of Analysis

No. CA06174

Product: Eichrom's Sr Resin Cartridges  
Lot #: SRSR6D

Eichrom's Sr Resin Cartridges, Lot # SRSR6D has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	124
Gravimetric Recovery	$\leq 105\%$	90.5 %
Chemical Recovery, Corrected	90-110 %	94.4 %



Prepared: December 20, 2006

Jill Bryant  
Quality Control Chemist