

Certificate of Analysis

No. CA07088

Product: Eichrom's Sr Resin Cartridges

Lot #: SRSR7B

Eichrom's Sr Resin cartridges, Lot # SRSR7B has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	133
Gravimetric Recovery	$\leq 105\%$	102 %
Chemical Recovery, Corrected	90-110 %	101 %



June 25, 2007

Jill Bryant
Quality Control Chemist