

Certificate of Analysis

No. CA08010

Product: Eichrom's Sr Resin cartridges
Lot #: SRSR8A

Eichrom's Sr Resin cartridge Lot # SRSR8A has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. The resin used to manufacture this lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	136
Gravimetric Recovery	$\leq 105\%$	104 %
Chemical Recovery, Corrected	90-110 %	94.7 %



January 23, 2008

Jill Bryant
Quality Control Chemist