

## Certificate of Analysis

No. CA08056

Product: Eichrom's Sr Resin cartridges


Lot #: SRSR8C

Eichrom's Sr Resin cartridge Lot # SRSR8C has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	116
Gravimetric Recovery	$\leq 105\%$	102 %
Chemical Recovery, Corrected	90-110 %	95.5 %

Jill Bryant  
Quality Control Chemist



March 25, 2008