

## Certificate of Analysis

No. CA08131

Product: Eichrom's Sr Resin 1mL cartridges  
Lot #: SRSR8E

Eichrom's Sr Resin 1mL cartridges, Lot # SRSR8E has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	111
Gravimetric Recovery	$\leq 105\%$	102 %
Chemical Recovery, Corrected	90-110 %	95.2 %



June 23, 2008

Jill Bryant  
Quality System Coordinator