

## Certificate of Analysis

No. CA05109

Product: Eichrom's TEVA Resin


Lot #: TEA10045

Eichrom's TEVA Resin, Lot # TEA10045, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

| <u>Test Value</u>                   | <u>Specification</u> | <u>Result</u> |
|-------------------------------------|----------------------|---------------|
| Appearance                          | White beads          | Pass          |
| Distribution Ratio, 3 N nitric acid | > 200                | 394           |

Signature



Prepared: October 20, 2005

Jill Bernstein

Quality Control Chemist