

Certificate of Analysis

No. CA05111

Product: Eichrom's TEVA Resin
Lot #: TEA10105

Eichrom's TEVA Resin, Lot # TEA10105, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3 N nitric acid	> 200	385

Signature Jill Bernstein
Jill Bernstein
Quality Control Chemist

Prepared: October 20, 2005