

expertise. commitment. results.

Certificate of Analysis

No. CA05112

Product: Eichrom's TEVA Resin

Lot #: TEA10135

Eichrom's TEVA Resin, Lot # TEA10135, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test ValueSpecificationResultAppearanceWhite beadsPassDistribution Ratio, 3 N nitric acid> 200371

Signature_

Prepared: October 28, 2005

Jill Bernstein Quality Control Chemist