

Certificate of Analysis

No. CA05116

Product: Eichrom's TEVA Resin


Lot #: TEA10275

Eichrom's TEVA Resin, Lot # TEA10275, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3 N nitric acid	> 200	440

Signature



Prepared: November 11, 2005

Jill Bernstein

Quality Control Chemist