

Certificate of Analysis

No. CA12129


Product: Eichrom's TEVA® Discs
Lot #: TEFD12A

Eichrom's TEVA® Discs, Lot # TEFD12A, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	47mm, round	Pass
Tc-99 yield in pH 2 (100-200 mL/min)	90 – 110 %	102.3 %

Jill Bryant
Quality System Coordinator



July 10, 2012