

Certificate of Analysis


No. CA08275

Product: Eichrom's TEVA® Discs
Lot #: TEFD8C

Eichrom's TEVA® Discs, Lot # TEFD8C, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	47mm, round	Pass
Tc-99 yield in pH 2 @100-200 mL/min	90 – 110 %	96.2 %



January 8, 2009

Jill Bryant
Quality Control Chemist