

**Certificate of Analysis**  
No. CA15091

Product: Eichrom's TEVA Resin  
Lot #: TES051215

Eichrom's TEVA Resin Lot # TES051215 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	474.0
% <sup>230</sup> Th stripped in 6M HCl		
Least acceptable mean (4 rep)	93.0	99.2
Max. standard deviation	2.5	0.1
% <sup>230</sup> Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.2
Max. standard deviation	1.5	0.1
% <sup>239</sup> Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	93.9
Max. standard deviation	2.5	1.7
% <sup>239</sup> Pu stripped in 6M HCl		
Maximum acceptable mean (4 rep)	1.0	0.2
Max. standard deviation	0.5	0.2
Extractant distribution ratio, <sup>230</sup> Th in 3M HNO <sub>3</sub>	>250	252
Extractant distribution ratio, <sup>230</sup> Th in 6M HCl	<0.01	0.001

  
\_\_\_\_\_  
Sarah Tejchma  
Quality Manager

*20 May 15*  
\_\_\_\_\_  
May 20, 2015