


Certificate of Analysis
 No. CA09161

Product: Eichrom's TEVA Resin
 Lot #: TES101909

Eichrom's TEVA Resin Lot # TES101909 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	339
% ²³⁰ Th stripped in 6M HCl		
Least acceptable mean (5 rep)	93.0	100.3
Max. standard deviation	2.5	1.6
% ²³⁰ Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (5 rep)	2.0	0.9
Max. standard deviation	1.5	0.3
% ²³⁹ Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (5 rep)	90.0	98.6
Max. standard deviation	2.5	0.3
% ²³⁹ Pu stripped in 6M HCl		
Maximum acceptable mean (5 rep)	1.0	0.4
Max. standard deviation	0.5	0.2



Joel Williamson
 Joel Williamson
 Quality Manager

Prepared: November 6, 2009