

Certificate of Analysis

No. CA06160

Product: Eichrom's TEVA Resin Cartridges


Lot #: TESR6A

Eichrom's TEVA Resin Cartridges, Lot # TESR6A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3 N nitric acid	> 200	398

Signature



Prepared: December 20, 2006

Jill Bryant
Quality Control Chemist