

**Certificate of Analysis**  
No. CA 05136


Product: Eichrom's TRU Resin  
Lot #: TRS10175

Eichrom's TRU Resin Lot # TRS10175 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

| Test Value                                                                     | Specification | Result |
|--------------------------------------------------------------------------------|---------------|--------|
| Appearance                                                                     | White beads   | Pass   |
| Distribution Ratio, Eu in 3.1 N nitric acid                                    | $\geq 60$     | 95     |
| % $^{241}\text{Am}$ stripped in 4M HCl                                         |               |        |
| Least acceptable mean (5 rep)                                                  | 95.0          | 99.5   |
| Max. standard deviation                                                        | 2.5           | 0.2    |
| % $^{241}\text{Am}$ stripped in 0.1M ammonium bioxalate                        |               |        |
| Maximum acceptable mean (5 rep)                                                | 1.0           | 0.2    |
| Max. standard deviation                                                        | 0.5           | 0.3    |
| % $^{239}\text{Pu}$ stripped in 0.1M ammonium bioxalate                        |               |        |
| Least acceptable mean (5 rep)                                                  | 95.0          | 96.8   |
| Max. standard deviation                                                        | 2.5           | 2.3    |
| % $^{239}\text{Pu}$ stripped in 4M HCl                                         |               |        |
| Maximum acceptable mean (5 rep)                                                | 1.0           | 0.2    |
| Max. standard deviation                                                        | 0.5           | 0.3    |
| Extractant distribution ratio, $^{241}\text{Am}$ strip in 0.05M $\text{HNO}_3$ | $<0.50$       | 0.24   |

TRU Resin Elution Curve



*Joel Williamson*

Joel Williamson  
Senior Quality Control Chemist

Prepared: December 22, 2005