

Certificate of Analysis

No. CA06014

Product: Eichrom's TRU Resin Cartridges

Lot #: TRSR5C

Eichrom's TRU Resin cartridge lot # TRSR5C has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 60	90.3



Jill Bernstein
Quality Control Chemist

Prepared: February 15, 2006