## SAFETY DATA SHEET

TRU Resin SDS - Bulk and Cartridge

Section 1: Chemical Product and Company Identification

Product Name TRU Resin

Product Number(s): TR1ML-R01-S, TR1ML-R50-S, TR-B01-A, TR-B01-F, TR-B01-S, TR-B01-XAD, TR-B100-A,

TR-B100-S, TR-B200-A, TR-B200-S, TR-B25-A, TR-B25-S, TR-B50-A, TR-B50-S, TR-R01-S,

Revision Date: 22-May-15

TR-R25-S, TR-R50-S

Product Synonym(s): Eichrom TRU Resin

Identified Uses: Laboratory chemicals, manufacture of substances

Manufacturer: Eichrom Technologies LLC General (8-5 CST M-F)

1955 University Lane Information: 800-422-6693 (in USA)

Lisle, Illinois 60532 630-963-0320

24 Hour Emergency Number:

CHEMTREC: 800-424-9300

### Section 2: Hazard(s) Identification



GHS Signal Word: Warning
GHS Classification of substance or mixture: Skin Irritant

Respiratory Tract Irritation

Acute toxicity, Inhalation (Category 5)

Acute hazards to the aquatic environment (Category 2)

Hazard Statement(s): H319 Causes serious eye irritation

H315 Causes skin irritation

H335 May cause respiratory irritation

H333 May be harmful if inhaled

H401 Toxic to aquatic life

| Pr         | -01 | 10  | nŧ | ia | n.  |
|------------|-----|-----|----|----|-----|
| $\Gamma$ I | E   | ٧U  | ш  | ıo | 11. |
|            | 0   | v C |    |    |     |

Disposal: P501

| Prevention:    |  |
|----------------|--|
| P261           | Avoid breathing dust and vapors.   |
| P264           | Wash hands thoroughly after handling.  |
| P271           | Use only outdoors or in a well-ventilated area.  |
| P273           | Avoid release to the environment.  |
| P280           | Wear protective gloves and eye protection.   |
| Response:      |  |
| P302/P352      | IF ON SKIN: Wash with plenty of soap and water.  |
| P304/P312      | Call a POISON CONTROL CENTER or doctor if you feel unwell.   |
| P304/P340      | IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing.  |
| P305/P351/P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P332/P313      | If skin irritation occurs, seek medical attention.   |
| P337/P313      | If eye irritation persists, get medical attention.   |
| P362           | Take off contaminated clothing and wash before reuse.  |
| P391           | Collect Spillage.  |
| Storage:       |  |
| P403/P233      | Store in a well-ventilated place. Keep container tightly closed.   |

Dispose of contents/container in accordance with federal, state, and local regulations.

Revision Date: 22-May-15

| Section 3: Composition / Information                                       |   |   |
|--|---|---|
| Component  | CAS_Number  | Percentage Range                            |
| Nonionic Acrylic Ester Polymer   |   | 60%   |
| TriButylPhosphate  | 126-73-8  | 27-28%                                      |
| CMPO   | 83242-95-9  | 12-13%                                      |
| Section 4: First-aid Measures  |   |   |
| Ingestion  | Contact local poison control center   |   |
| Skin Contact   | Wash immediately with soap and copious contaminated clothing promptly. If irritatio   |   |
| Eye Contact  | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.        |   |
| Inhalation   | Remove to fresh air. If breathing is labored, administer oxygen. If not breathing, give artificial respiration. Seek medical attention. |   |
| Most important symptoms and effects, both acute and delayed                | No further relevant information available.  |   |
| Indication of any immediate medical attention and special treatment needed | Treat according to symptoms (decontaminantidote.  | nation, vital functions), no known specific |
| Section 5: Firefighting Measures   |   |   |
| Extinguishing Media  | Foam, CO2, Dry Chemical   |   |
| Special Hazards  | Possible combustion products include phodioxide, and carbon monoxiede, in addition  |   |
| Fire and Explosion Hazards   | Highly toxic and irritating fumes may be rebe toxic.  |   |
|  | Polymer does not support flame.   |   |
| Protective Equipment   | Wear positive pressure self-contained bre-<br>protective equipment.   | athing apparatus and full personal          |
| Section 6: Accidental Release M  | easures   |   |
|  | Use proper personal protect equipment (s  | pecified in section 8)                      |
| Personal precautions   | Avoid breathing vapors, mist, or gas. See   |   |
|  | Surface may be slippery.  |   |
| Containment Cleanup  | Sweep up material and transfer to a suitab  | ole container for disposal.                 |
| Reference to other sections  | For disposal see section 13.  |   |
| Section 7: Handling and Storage  | •   |   |
| Specific End Use(s)  | Apart from the uses mentioned in section  | 1 no other specific uses are stinulated     |
| Conditions for safe handling   | Avoid contact with skin and eyes. Avoid in  |   |
| conditions for safe flanding   | Use mechanical exhaust if dust is formed.   | ·   |
| Conditions for safe storage  | Keep away from strong oxidizers.  |   |
| o en amono non o ano etenage   | Normal warehouse storage in cool, dry are   | ea is satisfactory                          |
| Section 8: Exposure Controls / F   |   | 54 15 54.16.46.6.y.                         |
| Engineering controls   | Mechanical exhaust is required.   |   |
| Control Parameters   | Per ACGIH, TLV-TWA for Tributylphospha  | ata in 2 nam                                |
| Control Farameters   | Per OSHA, PEL TWA is 5 mg/m3 for Tribu<br>Per NIOSH, REL TWA is 0.2 ppm (2.5 mg<br>ppm for Tributylphosphate                            | utylphosphate.                              |
| Exposure Controls  | Avoid contact with skin, eyes, and clothing immediately after handling the product.   | g. Wash hands before breaks and             |
|  | Do not eat, drink or smoke when using this  | s product                                   |
| Eye protection   | Wear safety glasses.  |   |
| Skin Protection  | Wear protective gloves, clothing, and eye   | protection.                                 |
| Respiratory protection   | Use NIOSH/MSHA approved respirator whechanical exhaust. An air-purifying respicanister may be permissible.                              |   |

Revision Date: 22-May-15

| Section 9: Physical Properties |   |   |                               |  |
|--------------------------------|---|---|-------------------------------|--|
| Information on basic ph        | Information on basic physical and chemical properties |   |                               |  |
| Appearance:                    | Powder  | Explosion Limits                          | Not Established               |  |
|                                | Hard, off-white spherical beads                       | (Upper/Lower):                            |                               |  |
| Odor:                          | Low to none   | Flash Point:                              | Not applicable                |  |
| Odor Threshold:                | Not Established                                       | Flammability:                             | Not Established               |  |
| pH:                            | Not relevant  | Autolgnition Temperature:                 | Product is not self-igniting. |  |
| Melting Point:                 | Not Established                                       | Decomposition Temperature Not Established |                               |  |
| Boiling Point:                 | Not Established                                       |   |                               |  |
| Relative Density:              | 0.37 g/mL at 25°C                                     | VaporPressure:                            | Not Established               |  |
| Solubility:                    | Insoluble in water                                    | VaporDensity:                             | Not Established               |  |
| Partition Coefficient:         | Not Established                                       | Evaporation Rate:                         | Not Established               |  |
| Viscosity:                     | Not Applicable  |   |                               |  |
| Section 10: Stability          | v and Reactivity                                      |   |                               |  |

| Section 10: Stability and Reactivity |  |  |
|--------------------------------------|--|--|
| Reactivity                           | No hazardous reactions if stored and handled as indicated.   |  |
| Chemical Stability                   | Stable under normal handling and storage conditions.   |  |
| Hazardous Reactions                  | No hazardous reactions are expected in normal laboratory use. Hazardous polymerization will not occur.   |  |
| Materials to Avoid                   | Contact with strong oxidizers will degrade material.   |  |
| Hazardous decomposition Products     | Possible combustion products include phosphorous oxides, phosphoric acid, carbon dioxide, and carbon monoxide; additional unidentified organic compounds may also be produced. |  |

# Section 11: Toxicology Information

| Acute Toxicity     |  |
|--------------------|--|
| Oral Effects       | CMPO, oral LD50 >5000 mg/kg.   |
|                    | TRU Resin, estimated oral LD50 is 4300 mg/kg (mouse).                        |
|                    | Polymer, Oral LD50 > 5,000 mg/kg (rat)                                       |
|                    | Tributylphosphate, Oral LD50 is 3000 mg/kg (rat).                            |
|                    | Tributylphosphate, Oral LD50 is 1189 mg/kg (mouse).                          |
| Inhalation Effects | Tributylphosphate, inhalation LC50, mouse: 1300 mg/m3.                       |
|                    | Tributylphosphate, inhalation LC50 is 2529 ppm (rat - 1 hr) [or 28000 mg/m3] |
|                    | TRU Resin, estimated inhalation LC50 is 9200 ppm (rat - 1 hr) [or 36 mg/l]   |
| Eye Effects        | Tributylphosphate, Draize test, rabbit, eye: 500 mg Severe.                  |
|                    | May cause irritation or corneal injury.                                      |
| Dermal Effects     | TRU Resin, estimated dermal LD50 is > 5,000 mg/kg (mouse).                   |
|                    | CMPO, dermal LD50 >2000 mg/kg.   |
|                    | Polymer, Dermal LD50 > 5,000 mg/kg (rabbit)                                  |
|                    | Tributylphosphate, dermal LD50 is >3100 mg/kg (rabbit)                       |

Skin corrosion/irritation

Tributylphosphate is a mild eye irritant (rabbit). No data available for other

components.

Serious eye damage/irritation

Tributylphosphate is a mild skin irritant (rabbit). No data available for other

components.

Respiratory or skin sensitization

No data available regarding respiratory or skin sensitization effects of this product.

Germ Cell Mutagenicity

No data available regarding mutagenic effects of this product.

Carcinogenicity

| TRU Resin SDS - Bulk and Cartridge | Revision Date: | 22-May-15 |
|------------------------------------|----------------|-----------|
|------------------------------------|----------------|-----------|

|                                | The result obe balk and saitings   |
|--------------------------------|--|
|                                | No specific data available for some components. Minimize direct exposure to material.  |
|                                | No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by IARC. |
|                                |  |
|                                | No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.                  |
|                                | No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.              |
|                                | Tributylphosphate, Limited evidence of carcinogenicity in animal studies. (Rat, Oral, Kidney, Ureter, Bladder: tumors; Mouse, Oral, Liver: tumors).      |
| Reproductive Toxicity          |  |
|                                | No data available regarding reproductive effects of this product.  |
| Specific Target Organ Toxicity |  |
| Single Exposure                | No data available regarding specific target organ toxicity single exposure.  |
| Repeated Exposure              | No data available regarding specific target organ toxicity repeated exposure.  |

No data available regarding aspiration hazards associated with this product.

#### Section 12: Ecological Information

Aspiration Hazard

| Cotton 12. Lociogica information        |   |  |  |
|---|---|--|--|
|   | No data are available on the adverse effects of this material on the environment.   |  |  |
| Aquatic Toxicity                        |   |  |  |
| Acute Toxicity to fish                  | Tributylphosphate - LC50, Carassius auratus (goldfish) - 8.8 mg/l - 96 hr   |  |  |
| Acute Toxicity to aquatic invertebrates | Tributylphosphate - EC50, Daphnia magna (water flea) - 3.6 mg/l - 48 hr   |  |  |
| Acute toxicity to aquatic plants        | Tributylphosphate - EC50, Desmodesmus subspicatus (green algae) - 1.1 mg/l - 72 hr  |  |  |
| Acute Toxicity to fish                  | TRU Resin - Estimated LC50, Carassius auratus (goldfish) - 32 mg/l - 96 hr  |  |  |
| Acute Toxicity to aquatic invertebrates | TRU Resin - Estimated EC50, Daphnia magna (water flea) - 14 mg/l - 48 hr  |  |  |
| Acute toxicity to aquatic plants        | TRU Resin - Estimated EC50, Desmodesmus subspicatus (green algae) - 4.5 mg/l - 72 hr  |  |  |
| Persistance and degradability           |   |  |  |
| Biodegradability                        | Tributylphosphate - aerobic biodegradability - exposure time 28 d, 89% - readily biodegradable.   |  |  |
| Bioaccumulative potential               |   |  |  |
|   | Tributylphosphate - Oryzias latipes - 38 d - 93 ug/l; Bioconcentration factor: 21-35  |  |  |
| Mobility in Soil                        | No data are available for mobility in soil.   |  |  |
| PBT/vPvB assessment                     | PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.   |  |  |
| Other                                   | An environmental hazard cannot be excluded in the event of unprofessional handling or disposal. Very toxic to aquatic life with long lasting effects. |  |  |

#### Section 13: Disposal Considerations

General Dispose of contents/container in accordance with federal, state, and local regulations.

Unused: Bury resin in licensed landfill or burn in approved incinerator equipped with an afterburner and scrubber according to local, state, and federal regulations.

Used: For resin contaminated with hazardous materials, dispose of mixture as hazardous material according to local, state, and federal regulations.

# Section 14: Transport Information

Air Transport: Not Hazardous per IATA 2014

Ground Transport: Not D.O.T. Hazardous

Water Transport: Not Hazardous per IMDG 2012.

#### Section 15: Regulatory Information

**US Federal Regulations** 

Toxic Substances Control Act (TSCA): This material is provided to you under the research and development (R&D) exemption.

**US State Regulations** 

A component, Tributylphosphate [CAS 126-73-8], is listed on the following state right to know lists:

CA, MA, MN, NJ, PA

Printed: Friday, May 22, 2015 Page 4 of 5 Eichrom Technologies LLC

# TRU Resin SDS - Bulk and Cartridge

## Section 16: Other Information

Revision Updated to GHS SDS format, including classification

SDS Prepared By: Eichrom Technologies LLC

The information set forth herein has been gathered from standard reference materials and is to the best knowledge and belief of Eichrom Technologies LLC, accurate and reliable. Such information is offered solely for your consideration, investigation and verification, and does not suggest or guarantee that the hazard precautions or procedures mentioned are the only ones that exist. Eichrom Technologies LLC makes no warrantees, express or implied, with respect to the use of such information or the use of the specific material identified herein in combination with any other material or process, and assumes no responsibility therefore.

Revision Date: 22-May-15