

Certificate of Analysis

No. CA06011

Product: Eichrom's UTEVA Resin

Lot #: UTA01196

Eichrom's UTEVA Resin, Lot # UTA01196, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 100	212

Signature Jill Bernstein
Jill Bernstein
Quality Control Chemist

Prepared: February 10, 2006