

## Certificate of Analysis

No. CA06036

Product: Eichrom's UTEVA Resin


Lot #: UTA03096

Eichrom's UTEVA Resin, Lot # UTA03096, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 100	227

Signature



Prepared: April 7, 2006

Jill Bryant  
Quality Control Chemist