

expertise. commitment. results.

Certificate of Analysis

No. CA06038

Product: Eichrom's UTEVA Resin

Lot #: UTS03166

Eichrom's UTEVA Resin, Lot # UTS03166, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	<u>Specification</u>	Result
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 100	226

Signature

Jill Bryant

Quality Control Chemist

Prepared: April 7, 2006