

Certificate of Analysis

No. CA06089

Product: Eichrom's UTEVA Resin

Lot #: UTS08156

Eichrom's UTEVA Resin, Lot # UTS08156, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 100	141

Signature Jill Bryant Prepared: August 22, 2006

Jill Bryant
Quality Control Chemist