

Certificate of Analysis


No. CA06193

Product: Eichrom's UTEVA Resin Cartridges
Lot #: UTSR6J

Eichrom's UTEVA Resin cartridges, Lot # UTSR6J, has been manufactured according to Eichrom standard operating procedures. The lot of resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 100	155



Jill Bryant
Quality Control Chemist

Prepared: January 5, 2007