



# Certificate of Analysis

No. CA15001

Product: Eichrom's Weak Base EC Resin

Lot #: WBS010615

Eichrom's WBEC® Resin Lot # WBS010615 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% <sup>233</sup> U Stripped in 1M NH <sub>4</sub> OH		
Maximum Acceptable Mean (4 Rep)	TBD	1.0
Maximum Standard Deviation	TBD	0.3
% <sup>99</sup> Tc Stripped in 1M NH <sub>4</sub> OH		
Least Acceptable Mean (4 Rep)	TBD	95.2
Maximum Standard Deviation	TBD	0.3
% <sup>233</sup> U stripped in 0.1M HNO <sub>3</sub>		
Least Acceptable Mean (4 Rep)	TBD	98.7
Maximum Standard Deviation	TBD	0.4
% <sup>99</sup> Tc stripped in 0.1M HNO <sub>3</sub>		
Maximum Acceptable Mean (4 Rep)	TBD	0.0
Maximum Standard Deviation	TBD	0.0

Sarah Tejchma  
QC Manager

16 Jan 15

January 16, 2015