

**New Chromatographic Materials for Determinations of Actinides,
Strontium, and Technetium in Environmental, Bioassay,
and Nuclear Waste Samples***

E. Philip Horwitz

Chemistry Division, Argonne National Laboratory, Argonne, IL

***Work performed under the auspices of the Office of Basic Energy Sciences, Division of
Chemical Sciences, U.S. Department of Energy, under contract number W-31-109-ENG-38.**

Co-workers

Chemical Separations Group, Chemistry Division

Mark Dietz

Renato Chiarizia (Visiting Scientist, ENEA)

Herbert Diamond

Cary Bauer (Student, Northern ILL Univ.)

Internal Dosimetry Group, Environment, Safety, and Health Division

Donald Nelson

Theresa Sullivan

Analytical Chemistry Group, Chemical Technology Division

Lesa Smith


TABLE OF CONTENTS

FUNDAMENTALS

A.	Schematic drawing of extraction chromatographic system	1
B.	Features of the new chromatographic materials.....	2
C.	Currently available extraction chromatographic materials.....	3
D.	Characteristics of extraction chromatographic materials and columns.....	4
E.	Relationship between solvent extraction (SX) and liquid-liquid chromatography (LLC)	5
F.	TEVA•Spec	
	1. Basic Chemistry.....	6
	2. k' of actinides III, IV, V, and VI vs. HNO_3 and HCl	7
	3. Effect of matrix constituents (anions) on Np(IV) retention	8
	4. Effect of matrix constituents (uranyl and thorium nitrates) on Np(IV) retention	9
	5. Separation of U(VI), Th(IV), and Np(IV).....	10
	6. Separation of Np and Pu.....	11
G.	U/TEVA•Spec	
	1. Basic Chemistry.....	12
	2. Characteristics of extraction chromatographic material and packed columns	13
	3. k' of actinides III, IV, and VI vs. HNO_3	14
	4. Elution behavior of selected elements.....	15
	5. Effect of matrix constituents (anions) on U(VI) retention	16
	6. Effect of matrix constituents (anions) on Np(IV) retention	17
H.	TRU•Spec	
	1. Basic Chemistry.....	18
	2. Characteristics of TRU•Spec extraction chromatographic material and packed columns	18a
	2. k' of actinides III, IV, V, and VI vs. HNO_3	19
	3. k' of selected ions vs. HNO_3	20
	4. Effect of matrix constituents on Am(III) retention from 2 M HNO_3	21
	5. Effect of Fe(II) and (III) on Am retention from 2 M HNO_3	22
	6. Effect of matrix constituents on Np(IV) retention from 2 M HNO_3	23
	7. Elution behavior of selected elements.....	24
	8. k' of actinides III, IV, and VI vs. HCl	25
	9. Effect of oxalic acid on U(VI) and Np(IV) retention from 1 M HCl	26
	10. Sequential separation of actinides	27
I.	Sr•Spec	
	1. Basic Chemistry.....	28
	2. Characteristics of extraction chromatographic material and packed columns.....	29
	3. k' of alkali and alkaline earths vs. HNO_3	30
	4. k' of actinides and selected ions vs. HNO_3	31
	5. Effect of matrix constituents on Sr retention from 3 M HNO_3	32
	6. Comparison of Sr elution curves for medium and small particle sizes.....	33
	7. Elution behavior of common elements and fission products on Sr•Spec.....	34

FUNDAMENTALS

EXTRACTION CHROMATOGRAPHY



$$k' = D V_s / V_m$$

k' = free column volumes to peak maximum

Features of the New Chromatographic Materials

- **High specificity**
- **Room temperature operation**
- **Gravity flow**
- **Acid - side loading from a wide range of concentrations**
- **Minimal waste generation**
- **Column compatibility for tandem arrangements**

Extraction Chromatographic Materials

	Currently Available	Specificity
Material		
TRU•Spec		Actinides (III, IV, VI), Ln(III)
U/TEVA•Spec		U(VI), Actinides (IV)
TEVA•Spec		Actinides (IV), Tc(VII)
Sr•Spec		Sr

Currently Under Development

Tc•Spec	Tc
Ra•Spec	Ra

**Characteristics of Extraction Chromatographic
Materials and Columns**

Stationary Phase **Neat Extractant or Extractant – Diluent**

Support **Amberlite XAD - 7**

Particle Sizes **50 - 100 μ .m, 100 - 125 μ .m**

Extractant Loading **40 weight percent**

Bed Density **0.33 to 0.37 g/mL**

Fractional Volumes:

Void **0.65**

Support **0.20**

Extractant **0.15**

Relationship Between SX and LLC

$$k' = D \cdot \frac{V_s}{V_m}$$

where

k' = retention volume

D = distribution ratio



V_s = volume of stationary phase

V_m = volume of mobile phase


Acid dependency of k' for various ions at 23°C.

TEVA • Spec


TEVA • Spec




Effect of Matrix Constituents on Neptunium Retention
TEVA • Spec/2 M HNO₃




Effect of Matrix Constituents on Neptunium Retention
TEVA • Spec/5 M HNO₃






Separation of U(VI), Th (IV), and Np (IV) on TEVA-Spec



Separation of Np and Pu on TEVA · Spec



U/TEVA·Spec

Extractant:**Diluent:** None**Extraction Equilibria:**

Characteristics of the U/TEVA-Spec extraction chromatographic material and packed columns.

	Bulk Material
Stationary Phase	Diamyl amyolphosphate ($\rho = 0.926$ g/ml)
Support	Amberchrom™ CG-71 or Amberlite™ XAD-7
Particle Diameter	50-100 μm (Amberchrom); 80-160 μm ; 100-125 μm
Extractant Loading	40 weight percent
Average Density of Extractant-Loaded Beads ^a	1.10 g/mL

Packed Columns

v_s	0.167 mL/mL of bed
Bed Density	0.386 g/mL
v_m (also FCV)	0.65 mL/mL of bed
v_s/v_m	0.257


Calculated Capacity

37 mg U/mL of bed

Experimentally Measured Capacity

37 mg U/mL of bed

^a Pycnometric density and flotation density were 1.065 and 1.133, using the 50-100 μm particle size resin. The calculated density is 1.094 g/ml assuming 100% pore filling and no swelling.

Acid dependency of k' for various ions at 23°C.

Elution behavior of selected elements on a U/TEVA-Spec (fine particle) column.

Element	Portion eluting (%) ^a						
	← 2 M HNO ₃ → 0.02 M HNO ₃ →						
	Number of free column volumes ^b						
	1-5	6-10	11-15	16-20	21-25	26-30	31-40
Li	99.4	—	—	—	—	—	—
Na	95.0	<1.2	—	—	—	—	—
Mg	101	—	—	—	—	—	—
Al	92.0	<5.8	—	—	—	—	—
K	<123	—	—	—	—	—	—
Ca	92.4	<13.3	—	—	—	—	—
Cr	94.7	<3.3	—	—	—	—	—
Mn	96.5	<0.7	—	—	—	—	—
Fe	94.0	<0.3	—	—	—	—	—
Co	99.1	<7.1	—	—	—	—	—
Ni	97.8	<2.0	—	—	—	—	—
Cu	98.1	<2.6	—	—	—	—	—
Zn	96.8	<1.2	—	—	—	—	—
Rb ^c	102	—	—	—	—	—	—
Sr	101	—	—	—	—	—	—
Y	97.9	<2.9	—	—	—	—	—
Zr	56.4	35.7	(1.4)	—	—	—	—
Ru	59.6	<19.2	<19.2	—	—	—	—
Rh	91.7	—	—	—	—	—	—
Ag	103	—	—	—	—	—	—
Cd	96.4	<0.7	—	—	—	—	—
Cs ^d	100	—	—	—	—	—	—
Ba	96.6	<14.4	—	—	—	—	—
La	108	—	—	—	—	—	—
Ce	96.3	<25	—	—	—	—	—
Pr	(136)	—	—	—	—	—	—
Nd	99.0	—	—	—	—	—	—
Sm	120	—	—	—	—	—	—
Eu	100	—	—	—	—	—	—
Pb	98.4	<5.8	—	—	—	—	—
U ^e	—	—	—	—	—	—	>99

^a Because of uncertainties inherent in the ICP-AES method used for quantitation, the fractions shown for a given element may not total 100%. Values in parentheses are subject to considerable uncertainty and are intended only as a guide. Feed solution contained ~ 0.02 M oxalic acid to solubilize zirconium.


^b 1 FCV = 0.60 mL here.

^c By flame atomic emission.


^d By atomic absorption.

^e Radiometric.

Effect of Matrix Constituents on Uranium Retention
UTEVA • Spec/2 M HNO₃




Effect of Matrix Constituents on Neptunium Retention
U/TEVA • Spec/2 M HNO_3





TRU·Spec

Extractant:



Diluent: TBP

Extraction Equilibria:



Characteristics of the TRU-Spec extraction chromatographic material and packed columns.

Bulk Material


Stationary Phase	0.75 M CMPO in TBP ($\rho = 0.971$ g/ml)
Support	Amberchrom™ CG-71
Particle Diameter	50-100 μm and 80-160 μm
Extractant Loading	40%
Average Density of Extractant-Loaded Beads ^a	1.12 g/mL

Packed Columns

v_s	0.152 mL/mL of bed
Bed Density	0.370 g/mL
v_m (also FCV)	0.68 mL/mL of bed
v_s/v_m	0.223
Calculated Capacity	5.49 mg Nd or 9.18 mg ²⁴¹ Am/mL of bed
Experimentally Measured Capacity	4.1 mg Nd or 6.8 mg ²⁴¹ Am/mL of bed


^a Pycnometric density and flotation density were 1.081 (in water) and 1.158 (in 4.9 M HNO₃) g/mL, respectively, using the 50-100 μm particle size resin. The calculated density is 1.094 g/ml assuming 100% pore filling and no swelling.

Acid dependency of k' for various ions at 23°C.




Acid dependency of k' for various ions at 23-25°C.

TRU • Spec




Effect of Matrix Constituents on Americium Retention TRU • Spec/2 M HNO₃



Effect of Fe on Am Uptake by TRU Spec
(Reference k'_{Am} at 2M HNO₃, 23°C = 98)

mg of Fe/10mL	M of Ascorbic Acid	k'_{Am}
25	-	28
50	-	17
100	-	9.1
25	0.3	88
50	0.3	72
100	0.3	65

Effect of Matrix Constituents on Neptunium Retention
TRU • Spec/2 M HNO_3



Elution of Selected Elements on a TRU-Spec Column
(*fine particles*)

Element	Fraction Eluting (%) by number of free column volumes						
	1-5	6-10	11-15	16-20	21-25	26-30	31-40
Li	98.4	<19	-	-	-	-	-
Na	92.8	<1.2	-	-	-	-	-
Mg	100	-	-	-	-	-	-
Al	99.8	<2.9	-	-	-	-	-
K	81.8	40.9	-	-	-	-	-
Ca	100	-	-	-	-	-	-
Cr	100	-	-	-	-	-	-
Mn	100	-	-	-	-	-	-
Fe	102	12.3	-	-	-	-	-
Co	100	-	-	-	-	-	-
Ni	100	-	-	-	-	-	-
Cu	100	-	-	-	-	-	-
Zn	100	-	-	-	-	-	-
Sr	100	-	-	-	-	-	-
Y	23.4	76.8	3.5	-	-	-	-
Zr	-	-	-	-	-	-	75.0
Ru	82.6	<19.2	-	-	-	-	-
Rh	100	-	-	-	-	-	-
Ag	100	-	-	-	-	-	-
Cd	100	-	-	-	-	-	-
Ba	100	-	-	-	-	-	-
La	-	-	-	-	-	30.0	72.0
Ce	-	-	-	-	-	<25.0	75.0
Pr	-	-	-	-	-	-	100
Nd	-	-	-	-	-	-	96.0
Sm	-	-	-	-	-	-	100
Eu	-	-	-	-	-	-	>99
Hg	[100]	[60]	[19]	-	-	-	-
Pb	100	-	-	-	-	-	-
Am*	-	-	-	-	-	-	>99


* radiometric

1-30 f.c.v. : 2M HNO₃


31-40 f.c.v. : 0.05M HNO₃


note: Because of uncertainties inherent in the ICP-AES method, the fractions shown for each element may not total to 100%. Values in parentheses are subject to considerable uncertainty and are intended only as a rough guide.

Acid dependency of k'
for various ions at 23°C.
TRU - Spec




TRU • Spec/HCl
Effects of oxalic acid on k'
of Np(IV) and U(VI) in 1 M HCl






Sr·Spec

Extractant: bis-t-Butyl-cis-Dicyclohexano-18-Crown-6



Diluent: 1-Octanol (C₈H₁₇OH)

Extraction Equilibrium:



CHARACTERISTICS OF STRONTIUM-SELECTIVE EXTRACTION
CHROMATOGRAPHIC MATERIAL AND PACKED COLUMNS

BULK MATERIAL		
Stationary Phase	1.0 M DtBuCH18C6 in 1-octanol ^a ($\rho = 0.912$ g/mL)	
Support	Amberchrom™ CG-71 or Amberlite™ XAD-7	
Particle Sizes (diam.)	50-100 μm ; 80-160 μm ; 100-125 μm	
Extractant Loading	40 weight percent	
Average Density of Extractant-Loaded Beads ^b	1.12 g/ml	
PACKED COLUMN		
	Particle Size	
	50-100 μm	80-160 μm
v_s , ml/ml of Bed ^c	0.146 \pm 0.004	0.152 \pm 0.004
Bed Density (g/ml)	0.33	0.35
v_m (Also F.C.V.), ml/ml of Bed ^c	0.71 \pm 0.01	0.69 \pm 0.01
v_s/v_m ^c	0.21 \pm 0.01	0.22 \pm 0.01
Capacity:		
Calc. mg Sr/ml of Bed	12.8	13.3
Exp. mg Sr/ml of Bed	10.6	—


^a 1 M concentration of DtBuCH18C6 in 1-octanol is close to a saturated solution; therefore, the Sr retention from nitric acid is close to the maximum that can be achieved with this system.

^b Pycnometric density and flotation density values were 1.079 and 1.157 g/ml, respectively, using the 50-100 μm particle size resin. The calculated density is 1.094 g/ml, assuming 100% pore filling and no swelling.

^c Average of five measurements: error is one σ .


Acid dependency of k' for various ions at 23-25°C.

Sr • Spec




Acid dependency of k' for various ions at 23-25°C.


Sr • Spec



Effect of Matrix Constituents on Strontium Retention
Sr • Spec/3 M HNO₃



Comparison of Elution Curves for Sr^{2+}
 for Two Particles of Sr · Spec
 Elutrient 3.2 M HNO_3 , 23-24°C



ELUTION BEHAVIOR OF COMMON ELEMENTS AND FISSION
PRODUCTS ON THE STRONTIUM-SELECTIVE RESIN

Element	PERCENT OF ELEMENT FOUND IN F.C.V. #						F.C.V.
	3 M HNO ₃ - 0.01 M Oxalic Acid						
	1-5	6-10	11-15	16-20	21-25	26-30	31-40
Li	100	—	—	—	—	—	—
Na	100	—	—	—	—	—	—
K	66	35	—	—	—	—	—
Rb	100	—	—	—	—	—	—
Cs	100	—	—	—	—	—	—
Mg	100	—	—	—	—	—	—
Ca	100	—	—	—	—	—	—
Sr	—	—	—	—	—	—	99
Ba	—	—	53	42	6	0.7	—
Ra	99		1	—	—	—	—
Al	100	—	—	—	—	—	—
Cr	100	—	—	—	—	—	—
Mn	100	—	—	—	—	—	—
Fe	99	0.6	0.2	0.4	—	—	—
Co	100	—	—	—	—	—	—
Ni	100	—	—	—	—	—	—
Cu	100	0.2	—	—	—	—	—
Zn	100	0.2	—	—	—	—	—
Y	100	0.1	—	—	—	—	—
Zr	91	0.4	0.2	—	—	—	—
Mo	84		16		—	—	—
Tc	57	43	—	—	—	—	—
Ru	100	—	—	—	—	—	—
Rh	100	—	—	—	—	—	—
Pd	100	—	—	—	—	—	—
Ag	15	88	2	—	—	—	—
Cd	100	0.1	—	—	—	—	—
La-Eu	100	0.1	—	—	—	—	—
Hg	5	5	19	40	19	10	5

Column parameters: Particle size = 50-100 μm , Bed Volume = 1.0 cm^3 , Bed height = 5.0 cm,
and 1 F.C.V. = 0.60 mL.