



# Certificate of Analysis

No. CA18021

Product: Eichrom's LN Resin  
Lot #: LNS022718  
Date of Manufacture: 27-Feb-2018

Eichrom's LN Resin Lot # LNS022718 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance Distribution Ratio, 0.1 N nitric acid	White beads >700	Pass 2589.7

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

March 1, 2018