



Certificate of Analysis

No. CA18016

Product: Eichrom's RE Resin
Lot #: REF021918
Date of Manufacture: 19-FEB-2018

Eichrom's RE Resin Lot # REF021918, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 125	154.0

Jill Bryant
Quality System Coordinator

February 22, 2018