

Certificate of Analysis No. CA20009

Product:

Eichrom's TRU Resin Cartridges

Lot #:

TRSR20A

Date of Manufacture: 16-JAN-2020

Eichrom's TRU Resin cartridge Lot # TRSR20A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, TRS010720, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	Result
Appearance	White beads	Pass
Distribution Ratio, Eu in 3.1 N nitric acid	≥ 60	101.9
% ²⁴¹ Am stripped in 4M HCl		
Least acceptable mean (4 rep)	95.0	99.5
Max. standard deviation	2.5	0.1
% 241Am stripped in 0.1M ammonium bioxala	ate	
Maximum acceptable mean (4 rep)	1.0	0.3
Max. standard deviation	0.5	0.1
% ²³⁹ Pu stripped in 0.1M ammonium bioxalat	e	
Least acceptable mean (4 rep)	95.0	97.9
Max. standard deviation	2.5	0.3
% ²³⁹ Pu stripped in 4M HCl		
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.0
Extractant distribution ratio, ²⁴¹ Am strip in 0.05M HNO ₃ < 0.50		0.17

Jill Bryant

Quality System Coordinator

January 17, 2020