

Certificate of Analysis

No. CA 21068

Product: Eichrom's Anion Exchange Resin Cartridges
Lot #: A8MR21B
Date of Manufacture: 5-MAR-2021

Eichrom's Anion Exchange resin cartridge Lot A8MR21B has been manufactured according to Eichrom standard operating procedures. This resin used to manufacture this lot, 031820A8M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
Wet Volume Capacity	≥ 1.2 meq/mL	1.4



Jill Bryant
Quality System Coordinator

March 5, 2021