

# Certificate of Analysis

No. CA21095

Product: Eichrom's TEVA® Discs  
Lot #: TEFD21A  
Date of Manufacture: 5-APR-2021

Eichrom's TEVA® Discs, Lot # TEFD21A has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	47mm, round	Pass
Tc-99 yield in pH 2 (100-200 mL/min)	90 – 110 %	94.4 %

  
Jill Bryant  
Quality System Coordinator

April 5, 2021