

Certificate of Analysis

No. CA 21161

Product: Eichrom's Anion Exchange Resin 1mL Cartridges

Lot #: A8MR21E

Date of Manufacture: 23-JUN-2021

Eichrom's Anion Exchange resin 1mL cartridge Lot A8MR21E has been manufactured according to Eichrom standard operating procedures. This resin used to manufacture this lot, 011921A8M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
Wet Volume Capacity	≥ 1.2 meq/mL	1.3



Sarah McAlister
Quality Manager

24 Jun 2021
June 24, 2021