

Certificate of Analysis

No. CA 21198

Product: Eichrom's Nickel Resin 5mL Columns
Lot #: NI21E
Date of Manufacture: 13-AUG-2021

Eichrom's Nickel Resin 5mL Columns Lot # NI21E has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, NIA052421, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	103.7 %
Ave. % Ni in Load/Rinse	<1%	0.2 %
Ave. % Ni in Strip	>99%	99.8 %



Jill Bryant
Quality System Coordinator

August 13, 2021