

Certificate of Analysis

No. CA 21205

Product: Eichrom's Sr Resin 5mL Columns

Lot #: SR21D

Date of Manufacture: 24-AUG-2021

Eichrom's Sr Resin 5mL Column Lot # SR21D has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, SRA072121, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value Appearance Distribution Ratio, Sr in 3.1 N nitric acid % Sr stripped in 0.05M HNO ₃	Specification White beads ≥ 100	Result Pass 119.1
Least acceptable mean (5 rep) Max. standard deviation	85.0 2.5	94.9 1.5
%Ca stripped in 0.05M HNO ₃ Maximum acceptable mean (5 rep) Max. standard deviation	0.25 0.25	0.03 0.01
%Ba stripped in 0.05M HNO ₃ Maximum acceptable mean (5 rep) Max. standard deviation	1.00 0.50	0.18 0.07
%Y stripped in 0.05M HNO₃ Maximum acceptable mean (5 rep) Max. standard deviation	0.25 0.25	0.0 0.0
%Y in 8M HNO₃ load and rinse Least acceptable mean (5 rep) Max. standard deviation	99.0 2.5	99.9 0.0

Jill Bryant

Quality System Coordinator

August 25, 2021

