

Certificate of Analysis

No. CA22058

Product: Eichrom's LN-3 Resin
Lot #: L3S030622
Date of Manufacture: 6-MAR-2022

Eichrom's LN-3 Resin Lot # L3S030622 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
La in Load/Rinse (%)	95-105	99.9
Be in Load/Rinse (%)	95-105	99.5
U in Load/Rinse (%)	< 3	0.2
La in HEDPA Strip (%)	< 1	0.1
Be in HEDPA Strip (%)	< 1	0.5
U in HEDPA Strip (%)	95-105	99.8



Jill Bryant
Quality System Coordinator

March 21, 2022