

Certificate of Analysis

No. CA22059

Product: Eichrom's LN-3 Resin Cartridges
Lot #: L3SR22A
Date of Manufacture: 18-MAR-2022

Eichrom's LN-3 Resin cartridge Lot L3SR22A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, # L3S030622, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
La in Load/Rinse (%)	95-105	99.9
Be in Load/Rinse (%)	95-105	99.5
U in Load/Rinse (%)	< 3	0.2
La in HEDPA Strip (%)	< 1	0.1
Be in HEDPA Strip (%)	< 1	0.5
U in HEDPA Strip (%)	95-105	99.8



Jill Bryant
Quality System Coordinator

March 21, 2022