

Certificate of Analysis

No. CA22043

Product: Eichrom's Tritium Columns
Lot #: H322A
Date of Manufacture: 18-FEB-2022

Eichrom's Tritium Columns, lot# H322A, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Radiometric yield of Tritium	90-110%	98.4



Jill Bryant
Quality System Coordinator

February 18, 2022