

Certificate of Analysis


No. CA22092

Product: Eichrom's LN Resin Columns
Lot #: LN22C
Date of Manufacture: 15-APR-2022

Eichrom's LN Resin column Lot LN22C has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, LNA052021, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2177



Jill Bryant
Quality System Coordinator

April 15, 2022