

Certificate of Analysis


No. CA 22075

Product: Eichrom's Sr Resin
Lot #: SRA032122
Date of Manufacture: 21-MAR-2022

Eichrom's Sr Resin Lot # SRA032122 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, Sr in 3.1 N nitric acid	≥ 100	120.5
% Sr stripped in 0.05M HNO ₃		
Least acceptable mean (5 rep)	85.0	90.1
Max. standard deviation	2.5	0.7
%Ca stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.00
Max. standard deviation	0.25	0.02
%Ba stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	1.00	0.12
Max. standard deviation	0.50	0.07
%Y stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.0
Max. standard deviation	0.25	0.0
%Y in 8M HNO ₃ load and rinse		
Least acceptable mean (5 rep)	99.0	100.0
Max. standard deviation	2.5	0.0


Jill Bryant
Quality System Coordinator

March 28, 2022