

## Certificate of Analysis No. CA22080

Product:

Eichrom's TEVA Resin

Lot #: Date of TES032922 29-Mar-2022

Manufacture:

Eichrom's TEVA Resin Lot # TES032922 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## **ANALYSIS RESULTS**

<u>Test Value</u>	Specification	Result
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid % <sup>230</sup> Th stripped in 6M HCl	>200	425.8
Least acceptable mean (4 rep)	93.0	99.4
Max. standard deviation	2.5	0.1
% <sup>230</sup> Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.0
Max. standard deviation	1.5	0.0
% <sup>239</sup> Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	96.8
Max. standard deviation	2.5	1.3
% <sup>239</sup> Pu stripped in 6M HCI		
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.1
Extractant distribution ratio, <sup>230</sup> Th in 3M HN0		267
Extractant distribution ratio, <sup>230</sup> Th in 6M HCl	<0.01	0.000

Jill Bryant

Quality System Coordinator

April 1, 2022

www.eichrom.com