

Certificate of Analysis No. CA22089

Product:

Eichrom's TEVA Resin Cartridges

Lot #:

TESR22F

Date of

12-APR-2022

Manufacture:

Eichrom's TEVA Resin cartridge Lot TESR22F has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, TES032922, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	<u>Specification</u>	Result
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	425.8
% ²³⁰ Th stripped in 6M HCI		
Least acceptable mean (4 rep)	93.0	99.4
Max. standard deviation	2.5	0.1
% ²³⁰ Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.0
Max. standard deviation	1.5	0.0
% ²³⁹ Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	96.8
Max. standard deviation	2.5	1.3
% ²³⁹ Pu stripped in 6M HCl	4.0	
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.1
Extractant distribution ratio, ²³⁰ Th in 3M HNO)3 >250	267
Extractant distribution ratio, ²³⁰ Th in 6M HCI	<0.01	0.000

Jill Bryant

Quality System Coordinator

April 13, 2022