

# Certificate of Analysis

No. CA 22114

Product: Eichrom's Pb Resin QML Cartridge  
Lot #: PBFR22C  
Date of Manufacture: 12-MAY-2022

Eichrom's Pb Resin QML Cartridge Lot PBFR22C has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, PBF070621, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.0 %
% Lead in Strip	> 99 %	100.0 %
Gravimetric Yield	80-120 %	102.6 %
Chemical Yield	80-120 %	88.5 %

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

May 16, 2022