

Certificate of Analysis

No. CA 22125

Product:

Eichrom's 25mm Resolve™ PE Filters in Disposable

Funnels

Lot #:

RFDF042522F

Date of Manufacture:

27-MAY-2022

Eichrom's 25mm Resolve™ PE Filters in Disposable Funnels Lot RFDF04252F has been manufactured according to Eichrom standard operating procedures. The PE filter membrane used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	<u>Specification</u>	Result
Am-241 FWHM avg keV (8 replicates)	<50.0	41.2
Am-241 FWHM values above 60 keV	0	0
Am-241 FWHM values between 50-60 keV	2 or fewer	0
Am-243 FWHM avg keV (8 replicates)	<50.0	39.3
Am-243 FWHM values above 60 keV	0	0
Am-243 FWHM values between 50-60 keV	2 or fewer	0

Jill Bryant

Quality System Coordinator

May 31, 2022

