

Certificate of Analysis

No. CA22131

Product:

Eichrom's TRU Resin

Lot #:

TRA060222

Date of Manufacture: 2-JUN-2022

Eichrom's TRU Resin Lot # TRA060222 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	Result
Appearance	White beads	Pass
Distribution Ratio, Eu in 3.1 N nitric acid	≥ 60	81.7
% ²⁴¹ Am stripped in 4M HCI		
Least acceptable mean (4 rep)	95.0	99.3
Max. standard deviation	2.5	0.3
% ²⁴¹ Am stripped in 0.1M ammonium biox	alate	
Maximum acceptable mean (4 rep)	1.0	0.4
Max. standard deviation	0.5	0.1
% ²³⁹ Pu stripped in 0.1M ammonium bioxa	late	
Least acceptable mean (4 rep)	95.0	98.2
Max. standard deviation	2.5	0.2
% ²³⁹ Pu stripped in 4M HCI		
Maximum acceptable mean (4 rep)	1.0	0.02
Max. standard deviation	0.5	0.1

Jill Bryant

Quality System Coordinator

June 9, 2022