

# Certificate of Analysis

No. CA 22173

Product: Eichrom's Nickel Resin Columns  
Lot #: NI22B  
Date of Manufacture: 21-JUL-2022

Eichrom's Nickel Resin column Lot NI22B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, NIA052522, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101.3 %
Ave. % Ni in Load/Rinse	<1%	0.2 %
Ave. % Ni in Strip	>99%	99.8 %

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

July 21, 2022