



Certificate of Analysis
No. CA 22139

Product: Eichrom's RE Resin
Lot #: REA061322
Date of Manufacture: 13-JUN-2022

Eichrom's RE Resin Lot # REA061322 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 125	148.7



Jill Bryant
Quality System Coordinator

June 15, 2022