

Certificate of Analysis



No. CA 22160

Product: Eichrom's Sr Resin 5mL Columns
Lot #: SR22C
Date of Manufacture: 14-JUN-2022

Eichrom's Sr Resin 5mL Column Lot # SR22C has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, SRA060122, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, Sr in 3.1 N nitric acid	≥ 100	112.9
% Sr stripped in 0.05M HNO ₃		
Least acceptable mean (5 rep)	85.0	97.8
Max. standard deviation	2.5	0.5
%Ca stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.00
Max. standard deviation	0.25	0.06
%Ba stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	1.00	0.34
Max. standard deviation	0.50	0.08
%Y stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.0
Max. standard deviation	0.25	0.0
%Y in 8M HNO ₃ load and rinse		
Least acceptable mean (5 rep)	99.0	99.9
Max. standard deviation	2.5	0.0


Sarah McAlister
Quality Manager
July 14, 2022