

Certificate of Analysis

No. CA 22184

Product: Eichrom's Anion Exchange Resin Cartridges
Lot #: A4MR22A
Date of Manufacture: 15-AUG-2022

Eichrom's Anion Exchange resin cartridge Lot A4MR22A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, # 091721A4M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d
Wet Volume Capacity	≥ 1.0 meq/mL	1.2


Jill Bryant
Quality System Coordinator

August 16, 2022