

Certificate of Analysis

No. CA 22179

Product: Eichrom's Actinide Resin 1mL Cartridges
Lot #: ACSR22E
Date of Manufacture: 8-AUG-2022

Eichrom's Actinide Resin 1mL cartridge Lot ACSR22E has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, ACS010721, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.5 M HCl acid	≥ 300	624.0



Jill Bryant
Quality System Coordinator

August 10, 2022