

Certificate of Analysis

No. CA 22190


Product: Eichrom's Pb Resin QML Cartridge
Lot #: PBFR22D
Date of Manufacture: 19-AUG-2022

Eichrom's Pb Resin QML Cartridge Lot PBFR22D has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, PBF070621, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.0 %
% Lead in Strip	> 99 %	100.0 %
Gravimetric Yield	80-120 %	102.6 %
Chemical Yield	80-120 %	88.5 %

Jill Bryant
Quality System Coordinator



August 19, 2022