

Certificate of Analysis

No. CA 22049

Product: Eichrom's Cation Exchange Resin; 50Wx8, 100-200 mesh, H+ Form
Lot #: 030122C8M
Date of Manufacture: 01-MAR-2022

Eichrom's Cation Exchange resin, Lot # 030122C8M has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Extractables	< 1 mg(residue) / g(resin)	0.12 mg/g
pH	≥ 4.0	pass
Wet Volume Capacity	min. 1.7 meq/mL	2.0 meq/mL (Per vendor cert)



Jill Bryant
Quality System Coordinator

March 3, 2022